

Low pH buffers Used to Separate Acidic Compounds in RP - HPLC Primer

Date: 22-APRIL-2012 Last Updated: 1-DECEMBER-2025

By suppressing their ionization at low pH, acids will be more like neutral compounds and will retain on Reversed Phase HPLC columns much the same way that neutral compounds will. The same is true for using bases with basic compounds.

Another benefit of running on Type A and Type B silica at low pH is that the silanol activity is reduced. This should sharpen your peaks.



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com