



## What are the dead volumes of the MRQ and Fused Insert autosampler vials - FAQ

For both MRQ and Fused Insert autosampler vials, the general residual volume is about 2  $\mu$ L.

*The dead volume is analogous to residual volume.*

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)