

## Are septa in autosampler vial caps self sealing inert and compatible with organic solvents - FAQ

**All of our non-slit septa are self-sealing and lined with Polypropylene, PTFE or other fluoropolymer for inertness.**

HPLC and GC solvents do not present a compatibility issue but acetonitrile and ethers are very volatile and may evaporate after a single use or piercing of the septa. You may need to replace the cap after an injection if you need extreme control over evaporation from the vial.

Click [HERE](#) for screw cap ordering information and pictures

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)