



## Have the amber autosampler vials changed color compared to earlier batches - FAQ

The amber vials catalog # *9502S-PP-AMBER* have been discontinued due to instabilities in the material which was causing color changes. These items have been replaced with catalog # [\*9502S-PP-AM\*](#). You may notice the newer vials are more opaque.

[Click HERE for 9502S-PP-AM Ordering Information](#)

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)