

Can BASIK brand screw caps be used in ion chromatography - FAQ

While BASIK™ brand screw caps are not designed for ion chromatography, they may work for some applications. But I would not use them sensitive methods.

We do suggest it might be worthwhile for maximum performance in IC to use our ultra-pure polyethylene screw caps, see image below.



Note that these caps, are a single use cap and will not reseal after injection.

Click <u>HERE</u> for more information about IC compatible caps. Click <u>HERE</u> for PP vials that are IC compatible.

Printed from the Chrom Resource Center Copyright 2025, All Rights Apply **MicroSolv Technology Corporation** 9158 Industrial Blvd. NE, Leland, NC 28451 Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com