



Using one syringe filter for many samples will I get sample carry over - FAQ

You will most likely get some kind of carry over. It is not recommended to use one filter for more than one sample.

[Click HERE for AQ syringe filter ordering information and pictures](#)

Attachment

MicroSolv filters equivalency study pdf [Download File](#)

Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com