


Is the MicroSolv AQ brand nylon syringe filter compatible with formamide - FAQ


Yes, our MicroSolv AQ™ nylon syringe filters are compatible with formamide.

Formamide, dimethylformamide, and acetamide, are all compatible with our MicroSolv AQ™ brand nylon syringe filters.


Acetamide:



Dimethylformamide:



Formamide:



Click [HERE](#) for ordering information and pictures of AQ brand of syringe filters.

Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com