

What is the correct thread type and tubing to use with the filter sparging stones - FAQ

Our filter / sparging stones have an internal thread of 5-40. These should be used with low pressure 1/8" or 1/16" ID polymeric tubing. The threads are actually small barbs that grab onto the tubing, ensuring a secure fit.

NOTE: It may be important to use the internal threaded sparging stones instead of those with fittings to avoid the mixture of many different components in the solvent. This type of connection of the tubing and the stone should not come off in situ as can happen with the tubing/stem configuration.

Click [HERE](#) for filter sparging stone ordering information and pictures

Click [HERE](#) for low pressure tubing ordering information and pictures

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com