



Which HPLC column ID dimensions are the Hichrom brand guard column holders cat no 82000-GH compatible with - FAQ

The Hichrom brand guard column holder with the catalog number [82000-GH](#) can be used with any guard column in the range of 1.0mm – 4.6mm internal diameter (ID). We do not suggest using these guard columns with analytical columns outside of this size range.

In addition to the holder and guard column, you will need a coupler which is used to connect the guard column holder to the analytical column. Click [HERE](#) for ordering information and a picture of the union / column coupler.

For smaller ID columns we recommend using a Cogent Column Filter instead such as 49650-005. Click [HERE](#) for ordering information for 2um Cogent Column Filters.



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com