



Does the direction of a guard column in the holder matter - FAQ

No, the guard column cartridge can be placed in the holder in either direction.

The label on the column does not indicate a particular direction to place the column in the holder. The guard columns have no flow direction in the holder.

[Click HERE for Ordering Information Including Guard Column Holders](#)



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com