PRODUCT SHEET



Slurry Reservoir Kit for 2.1 mm ID columns 30 mm & shorter 1. Slurry Reservoir, Black PEEK 4.6 mm ID, 50 mm long 2. End fitting to connect tubing from a pump to reservoir 3. PEEK sealing rings for the bottom and top connection 4. Packing adapter to connect the reservoir to the column 5. ETFE sealing ring for the top of the column connection

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451 tel. (732) 380-8900, fax (910) 769-9435

Email: info@mtc-usa.com Website: www.mtc-usa.com